


# Solid hydrogen micro spheres for laser proton acceleration



Jens Polz

Georg Becker, Rui Costa Fraga, Robert Grisenti, Anton Kalinin, Diethard Klöpffel,  
Alex Robinson and Malte Kaluza

Institut für Optik und Quantenelektronik, Friedrich-Schiller-Universität, Jena, Germany  
Helmholtz Institut Jena, Friedrich-Schiller-Universität, Jena, Germany  
Institut für Kernphysik, Goethe-Universität, Frankfurt a. M., Germany  
Plasma Physics Group, Central Laser Facility, STFC Rutherford-Appleton Lab., UK



Pellets too large for laser proton acceleration  
 $\sim 50 \mu\text{m}$  diameter

Pellet target facility too large for easy  
 implementation  $\sim 4.5 \text{ m}$

Low spatial ( $\sim \text{mm}$ ) and temporal stability  
 due to travel distance of pellets ( $\sim 2.5 \text{ m}$ )

**Not suited for laser plasma experiments**

**→ Own development**



## Space!

limiting size of cryostat

- limiting cooling power at operating temperature ( $\sim 5$  W)
- limiting gas flow rates of  $H_2$
- limiting nozzle diameter ( $5 \mu\text{m}$ )



droplet generator limited to size  $\sim 20$  cm due to spatial constraints in the experimental chamber

## Stability


synchronizable to the laser

temporal and spatial stability of the order of the laser focal spot diameter ( $\sim 3 \mu\text{m}$ )

When injected into vacuum, strong evaporation of H<sub>2</sub> and Ar leads to rapid cooling and freezing of jets



Additional volume with low pressure needed to slow down freezing




6 cm





Additional volume with low pressure needed to slow down freezing

Second capillary with larger pinhole (150  $\mu\text{m}$ ) providing suitable atmosphere to prevent freezing





Additional volume with low pressure needed to slow down freezing

Second capillary with larger pinhole ( $150 \mu\text{m}$ ) providing suitable atmosphere to prevent freezing




Piezo actuator triggering even breakup of jet to droplets












**Argon**





**Hydrogen**










- Transparent stuff is evident
- Two component Oxford glue might flow up to the capillary inlet. Severe issue especially if the entrance surface not flat (increases the liquefied-gas/glue contact surface)



## Target requirements for generation of mono energetic protons

- limited proton source size






## Target requirements for generation of mono energetic protons

- second ion species




## Target requirements for generation of mono energetic protons


- second ion species



shoot mass limited, multi species targets (e.g.  $H_2/D_2$  mixtures)  
and get mono energetic protons for free!



# Summary



**Compact source for solid hydrogen micro spheres suitable for laser plasma experiments**

**Experiment scheduled for next month**

**Exploration of benefits from mass limited pure hydrogen targets**

**Controlled adding of heavier ion species ( $D_2$ ) for mono energetic proton generation**